

Rocket[®] Craft[™] Suction Pumps

The Rocket® Craft™ Suction Pump and the Rocket® DUO-VAC™ Suction Pump have been developed to provide smooth, low volume vacuum at pre-determined negative pressures, so allowing simple, low turbulence, oocyte recovery

- HIGHLY RELIABLE: Rocket® Craft™ Pumps have a proven 20 year history of reliability and excellent value.
- NEW SINGLE USE WATER TRAP SET provides cost effective solution to regular processing of reusable water traps and tubing.
- SIMPLY, SAFE OPERATION with air-operated foot controls. High vacuum foot control
 available on Duo-Vac Suction Pumps.
- SUPER QUIET diaphragm pumps make it easy to live with in the OR and means low maintenance service.
- FULLY VARIABLE VACUUM from 0-400mmHg.



Ordering Information		
Rocket® Craft™ Suction Pump	240V/50Hz	R29655
Rocket® Craft™ Suction Pump	120V/60Hz	R29654
Rocket® Craft™ DUO-VAC Suction Pump	240V/50Hz	R29660
Rocket® Craft™ DUO-VAC Suction Pump	120V/60Hz	R29661
Rocket® Water Trap Set suitable for all Craft™ Suction Pumps packed individually, 10 per carton	Sterile each	R57.685

